

Dragon – What Can I Say?



Common Commands

New Lines and Paragraphs

1. New Line – starts a new line
2. New Paragraph – starts a new paragraph

Formatting Commands

1. All Caps [word]
2. Bold [word]
3. Cap [word]

Editing Commands

1. Delete [word] (or) Scratch [word]
2. Delete That (or) Scratch That – deletes the last spoken phrase
3. Insert Before [word]
4. Insert After [word]
5. Undo That – Executes cntl + Z in the application
6. Redo That – Executes cntl + Y in the application
7. Resume with [word] – deletes text up to the given word

Navigating Commands

1. Go to Bottom – moves the cursor to the bottom of the field
2. Go to End of Paragraph – moves the cursor to the end of the paragraph
3. Go to End of Sentence – moves the cursor to the end of the sentence
4. Go to Top – moves the cursor to the top of the field
5. End of Sentence

Selecting Text:

1. Select [word] – select a word
2. Select [word] through [word] – select multiple words
3. Select and Replace – select text and dictate over selected text
4. Go Back – returns the cursor to the previous point after the Select and Replace has been executed
5. Unselect That – unselect a selected word



Controlling Dragon

Tell Dragon “What Can I Say” and Dragon will display a list of available voice commands grouped by functionality and describes what you action each command executes.

1. Manage Auto-Texts – Opens the Manage Auto-Texts dialog box
2. Manage Commands – Opens the Manage Commands dialog box
3. Manage Vocabulary – Displays the Manage Vocabulary page
4. Manage Words – Displays the Manage Words page
5. Open Dragon – Maximizes the application and dictation box
6. Close Dragon – Minimizes the application to the task bar
7. Open Dictation Box – Opens the dictation box
8. Close Dictation Box – closes the dictation box
9. Discard Text – Deletes text in the dictation box
10. Launch Help – Opens the Dragon Help
11. Quit Help – Closes Dragon Help
12. Recall Text – Re-writes text in the dictation box
13. Transfer Text – Transfers text to target application

